
Uricase rekombinant E. coli, min. 4 U/mg
lyophil., E.C. 1.7.3.3

Cat.No. : 37830

Parameter	Method	Specification
Molecular weight		ca. 125 000
Appearance		light brown lyophilisate
Activity (U/mg)		min. 4
Storage (°C)		-15 to -25

Unit definition

1 unit is that amount of enzymatic activity which catalyzes the conversion of 1 μ mole uric acid per minute at 25°C, pH 8.5 to allantoin.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.