

Uridin-5'-triphosphat·Na₃-salz reinst	Kat.Nr./Cat.No. :	37960
Uridine-5'-triphosphate·Na ₃ -salt research grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₉ H ₁₂ N ₂ O ₁₅ P ₃ ·Na ₃ ·2H ₂ O
Molekulargewicht molecular weight		586.2
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	HPLC	min. 91.0 (UTP)
Trocknungsverlust (%) loss on drying	5 h, 105 °C	max. 8.0
Schwermetalle (ppm) heavy metals	als /as Pb	max. 10
IR		entspricht corresponds
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.