
Valinomycin reinst Valinomycin research grade	Kat.Nr./Cat.No. :	38076
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₅₄ H ₉₀ N ₆ O ₁₈
Molekulargewicht molecular weight		1111.3
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	HPLC	93.0 – 100.0
[α]_D/20°C	1.6% in Benzol 1.6% in benzene	+30° bis/ to +32°
Identität identity	IR	entspricht corresponds
Schmelzpunkt (°C) melting point		187 – 191
Wasser (%) water	KF	max. 0.5
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.