

---

<b>Wortmannin aus <i>Penicillium fusiculosum</i> reinst</b> Wortmannin from <i>Penicillium fusiculosum</i> research grade	<b>Kat.Nr./Cat.No. :</b>	<b>38385</b>
--	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>23</sub> H <sub>24</sub> O <sub>8</sub>
<b>Molekulargewicht</b> molecular weight		428.4
<b>Aussehen</b> appearance		weißes Pulver white powder
<b>Löslichkeit</b> solubility	5mg/ml in Methanol	klare farblose Lösung clear colorless solution
<b>Gehalt (%)</b> assay	TLC, HPLC	min. 97.5
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.