

Xylencyanol FF **Kat.Nr./Cat.No. : 38505**
C.I. 43535

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{25}H_{27}N_2O_7S_2 \cdot Na$
Molekulargewicht molecular weight		554.60
Aussehen appearance		dunkelrotes bis dunkelblaues Pulver dark red to dark blue powder
λ max. (nm)	0.001% in H ₂ O	613 - 617
A 1cm/λmax.	0.001% in H ₂ O	min. 0.8
Wasser water	KF	max. 10.0
TLC	Kieselgel 60 F ₂₅₄ n-Butanol 5 Ethanol 2 Eisessig/ glacial acetic acid 1 Wasser/ water 2	einheitlich, Spuren möglich one major spot, traces possible
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.