

Zeatin ribosid reinst Zeatin riboside research grade trans isomer	Kat.Nr./Cat.No. :	38568
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₅ H ₂₁ N ₅ O ₅
Molekulargewicht molecular weight		351.36
Aussehen appearance		weißes bis fast weißes krist. Pulver white to off-white cryst. powder
Gehalt (%) assay	TLC	min. 96.0
Wasser (%) water	KF	max. 1.0
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.