

---

**Protease-Inhibitor-Mix B**

**Cat.No. : 39105**

---

<b>Cat. number</b>	<b>Components</b>	<b>Target Protease</b>
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
11280	EDTA-Disodium	Metalloproteases
<b>Storage (°C)</b>		-15 to -25
<b>Reconstitution</b>	Dissolve the content of the vial (26.4mg) in 1 ml 90% DMSO to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). Slight turbidity in the reconstituted solution is normal. After reconstitution, aliquot and freeze at -20°C (stable for at least 4 weeks).	

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.