
Protein Test Mixture 4

Cat.No. : 39208

Protein components	Molecular weight
---------------------------	-------------------------

Carbonic anhydrase	29 000
---------------------------	--------

Albumin egg	45 000
--------------------	--------

Albumin bovine	67 000
-----------------------	--------

Phosphorylase B	97 400
------------------------	--------

Storage (°C)	+2 to +8
---------------------	----------

Reconstitution

Dissolve to a concentration of 1mg/ml in SDS-sample buffer.
A volume of 5µl is sufficient for SDS-PAGE (in a standard
form of 80 x 60 x 0.75mm).

SDS-sample buffer:

125 mM Tris-HCl pH 6.8; 2% SDS; 15% glycerol;
10 mM DTT; 0.025% Bromophenol blue; 0.025% Orange G

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.
--