

SERVA Unstained SDS PAGE Protein Marker
6.5 – 200 kDa, Liquid Mix

Cat.No. : 39215

Protein Components	Molecular weight (kDa)
Myosin (bovine muscle)	200
β-Galactosidase (E. coli. recombinant)	116
Albumin (bovine)	67
Albumin (egg)	45
Carbonic anhydrase (bovine erythrocytes)	29
Trypsin inhibitor (soybean)	21
Lysozym (egg)	14.3
Aprotinin (bovine lung)	6.5
Storage (°C)	-15 to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.