

SERVA Recombinant SDS PAGE Protein Marker
10 – 150 kDa, Liquid Mix

Cat.No. : 39218

Protein recombinant	Molecular weight (kDa)
Protein 1	150
Protein 2	100
Protein 3	80
Protein 4	60
Protein 5	40
Carbonic anhydrase, covalently stained	approx. 30
Protein 6	20
Protein 7	10
Storage (°C)	-15 to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.