

n-Butanol molecular biology grade	Kat.Nr./Cat.No. :	39550
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₄ H ₁₀ O
Molekulargewicht molecular weight		74.12
Aussehen appearance		klare farblose Flüssigkeit clear colorless liquid
Gehalt (%) assay	GC	min. 99.5
Wasser (%) water	KF	max. 0.1
2-Butanol (%)		max. 0.1
Dibutylether (%)		max. 0.1
Ketone (%) ketones		max. 0.02
Brechungsindex refractive index	n _D ²⁰	1.3995 ± 0.0005
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.