

Calciumchlorid · 2H₂O	Kat.Nr./Cat.No. :	39551
Calcium chloride · 2H ₂ O molecular biology grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		CaCl ₂ ·2H ₂ O
Molekulargewicht molecular weight		147.0
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	Titration	min. 99.0
pH	5% in H ₂ O	4.5 – 9.2
Sulfat (%) sulfate	SO ₄	max. 0.03
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.001
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.