

Ethanol unvergällt absolut Ethanol undenaturated absolute molecular biology grade	Kat.Nr./Cat.No. :	39556
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂ H ₆ O
Molekulargewicht molecular weight		46.07
Aussehen appearance		klare farblose Flüssigkeit clear colorless liquid
Gehalt (%) assay	GC	min. 99.7
Wasser (%) water	KF	max. 0.2
Brechungsindex refractive index	n _D /20°C	1.360 – 1.362
Identität identity	IR	entspricht corresponds
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.