

Isoamyl alcohol molecular biology grade	Kat.Nr./Cat.No. :	39557
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₅ H ₁₂ O
Molekulargewicht molecular weight		88.15
Aussehen appearance		klare farblose Flüssigkeit clear colorless liquid
Gehalt (%) assay	alle Isomere total isomers	min. 99.0
Säuregehalt (%) acidity	as CH ₃ COOH	max. 0.003
Farbe / color	APHA	max. 10
Schwermetalle (%) heavy metals	als/ asPb	max. 0.00001
Wasser (%) / water	K.F.	max. 0.5
Zink (%) / zinc	Zn	max. 0.00001
Abdampfrückstand (%) residue on evaporation		max. 0.002
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Protease Aktivität protease activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.