
Isopropanol molecular biology grade	Kat.Nr./Cat.No. :	39559
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₃ H ₈ O
Molekulargewicht molecular weight		60.09
Aussehen appearance		klare farblose Flüssigkeit clear colorless liquid
Gehalt (%) assay	GC	min. 99.7
Dichte (g/ml) specific density	20°C	0.785 – 0.786
Brechungsindex refractive index	20°C	1.376 – 1.378
Wasser (%) water		max. 0.1
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.