

Kaliumacetat molecular biology grade Potassium acetate	Kat.Nr./Cat.No. : 39567
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂ H ₃ O ₂ K
Molekulargewicht molecular weight		98.14
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration	min. 99.0
pH	5% in H ₂ O	7.5 – 9.0
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.001
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.