
Denhardt's Solution, lyophilisiertes Pulver Denhardt's solution, lyophilized powder (for 50x concentrate) molecular biology grade	Kat.Nr./Cat.No. :	39602
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		cremeweißes kristallines Pulver creamwhite crystalline powder
Löslichkeit solubility	30g/l in H ₂ O	klare farblose Lösung clear colorless solution
pH		6.8 – 7.8
Zusammensetzung composition		
Albumin Fraktion V		0.5g/ 50ml
Polyvinylpyrrolidone		0.5g/ 50ml
Ficoll® 400		0.5g/ 50ml
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.