

Formamid Formamide molecular biology grade	Kat.Nr./Cat.No. :	39604
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		CH ₃ NO
Molekulargewicht molecular weight		45.04
Aussehen appearance		farblose Flüssigkeit colorless fluid
Gehalt (%) assay	aus/ from N	min. 99.5
Identität identity	IR	entspricht corresponds
Wasser (%) water	KF	max. 0.1
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.0001
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.