

---

<b>Polyvinylpyrrolidon 90</b> molecular biology grade	<b>Kat.Nr./Cat.No. :</b>	<b>39609</b>
--	--------------------------	--------------

---

Parameter parameter	Methode method	Spezifikation specification
<b>Gewichtsmittelwert</b> average of weight		ca. 1 100 000
<b>Aussehen</b> appearance		weißes Pulver white powder
<b>Feststoffgehalt (%)</b> solid content		min. 95
<b>K-Wert</b> K-value		88 – 96
<b>Sulfatasche (%)</b> sulfated ash		max. 0.02
<b>Vinylpyrrolidon (ppm)</b> vinylpyrrolidone	HPLC	max. 100
<b>DNase/ RNase Aktivität</b> DNase/ RNase activity		nicht nachweisbar none detected
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> We do not guarantee that the product can be used for a special application.</p>
--