

Ammoniumacetat molecular biology grade Ammonium acetate	Kat.Nr./Cat.No. : 39750
---	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂ H ₇ NO ₂
Molekulargewicht molecular weight		77.08
Aussehen appearance		farbloses kristallines Pulver colorless crystalline powder
Gehalt (%) assay	Titration	min. 97.0
pH	5% in H ₂ O	6.5 – 7.3
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.001
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.