

Ammoniumacetat-Lösung 7.5M Ammonium acetate solution 7.5M molecular biology grade	Kat.Nr./Cat.No. :	39751
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		klare farblose Lösung clear colorless solution
Zusammensetzung composition		
Ammoniumacetat (#39750) ammonium acetate		578.1 g/l
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.001
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.