

Ammoniumchlorid Ammonium chloride molecular biology grade	Kat.Nr./Cat.No. : 39752
---	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		NH ₄ Cl
Molekulargewicht molecular weight		53.5
Aussehen appearance		farbloses krist. Pulver colorless cryst. powders
Gehalt (%) assay	Titration	min. 99.5
pH	5% in H ₂ O	4.5 – 5.5
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.0005
Eisen (%) lead	Fe	max. 0.0002
Sulfat (%) sulfate	SO ₄	max. 0.002
Phosphat (%) phosphate	PO ₄	max. 0.0002
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.