
Dithioerythrit DTE molecular biology grade	Kat.Nr./Cat.No. :	39758
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₄ H ₁₀ O ₂ S ₂
Molekulargewicht molecular weight		154.3
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration GC	min. 99.0
Schmelzpunkt (°C) melting point		82 - 86
A 1cm/ 283nm	0.02M in H ₂ O	max. 0.04
Identität identity	IR	entspricht corresponds
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.