

Ethylenglycol-bis(2-aminoethylether)-N,N,N',N'-tetraessigsäure (EGTA) Ethylene glycol bis(2-aminoethylether)-N,N,N',N'-tetraacetic acid molecular biology grade	Kat.Nr./Cat.No. :	39762
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₄ H ₂₄ N ₂ O ₁₀
Molekulargewicht molecular weight		380.35
Aussehen appearance		feines, weißes Pulver fine, white powder
Gehalt (%) assay	Titration	min. 99.0
Identität identity	IR	entspricht corresponds
Trocknungsverlust (%) loss on drying		max. 2.0
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.