

Kaliumchlorid Potassium chloride molecular biology grade	Kat.Nr./Cat.No.:	39768
--	-------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		KCl
Molekulargewicht molecular weight		74.6
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	Titration	min. 99.0
Trocknungsverlust (%) loss on drying		max. 1.0
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 5
Eisen (ppm) iron		max. 10
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.