

| | | |
|---|--------------------------|--------------|
| Magnesiumchlorid-6H₂O Magnesium chloride·6H ₂ O molecular biology grade | Kat.Nr./Cat.No. : | 39771 |
|---|--------------------------|--------------|

| Parameter parameter | Methode method | Spezifikation specification |
|--|--------------------------|--|
| Summenformel formula | | MgCl ₂ ·6H ₂ O |
| Molekulargewicht molecular weight | | 203.3 |
| Aussehen appearance | | farblose Kristalle colorless crystals |
| Gehalt (%) assay | Titration | 98.0 – 101.0 |
| Löslichkeit solubility | 10 % in H ₂ O | klar und farblos clear and colorless |
| Sulfat (%) sulfate | SO ₄ | max. 0.01 |
| Calcium (%) | Ca | max. 0.1 |
| Schwermetalle (%) heavy metals | als/ as Pb | max. 0.001 |
| DNase/ RNase Aktivität DNase/ RNase activity | | nicht nachweisbar none detected |
| Lagerung (°C) storage | | +15 bis / to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.