

Magnesiumsulfat heptahydrat Magnesium sulfate heptahydrate molecular biology grade	Kat.Nr./Cat.No. :	39773
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		MgSO ₄ · 7H ₂ O
Molekulargewicht molecular weight		246.48
Aussehen appearance		weißes kristallines Pulver white crystalline powder
Gehalt (%) assay	Titration	min. 99.5
pH	5% in H ₂ O	5.0 – 9.2
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.001
Chlorid (%) chloride	Cl	max. 0.0005
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.