

Piperazin-N,N'-bis(2-ethansulfonsäure) (PIPES) Piperazine-N,N'-bis(2-ethane sulfonic acid) molecular biology grade	Kat.Nr./Cat.No. :	39779
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_8H_{18}N_2O_6S_2$
Molekulargewicht molecular weight		302.4
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titr.	min. 99.0
A 1cm/ 250 nm	2% in NaOH	max. 0.2
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.0005
Identität identity	IR	entspricht corresponds
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.