

di-Natriumhydrogenphosphat·2H₂O di-Sodium hydrogen phosphate·2H ₂ O molecular biology grade	Kat.Nr./Cat.No. :	39783
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		Na ₂ HPO ₄ ·2H ₂ O
Molekulargewicht molecular weight		178.0
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	Titration	min. 99.5
pH	5% in H ₂ O	9.0 - 9.2
Gesamtstickstoff (%) total nitrogen	N	max. 0.001
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.001
Ursprung origin		nicht tierisch non animal
DNase/ RNase Aktivität DNase/ RNase activity		nicht nachweisbar none detected
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.