

SPEZIFIKATION
SPECIFICATION

SERDOLIT® CG-120 I p.A.
100 – 200 mesh analytical grade

Kat.Nr./Cat.No.:

40311

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		beiges Korn beige particles
Gegenion counter ion		Na ⁺
Trocknungsverlust (%) loss on drying	3 h, 105°C	50.0 – 60.0
Kapazität (eq/l) capacity		min. 1.3
Lagerung (°C) storage		+15 bis/ to +30

SERVA Electrophoresis GmbH