

Serdolit® MB-1 p.A. Serdolit® MB-1 analytical grade	Kat.Nr./Cat.No. :	40701
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		braunes transparentes Korn brown-colored transparent particle
Zusammensetzung composition		Mischbettionenaustauscher aus stark saurem KIA, H ⁺ und stark basischem AIA, Typ I, OH ⁻ , 1:1,5 mixed-bed ion exchanger composed strongly acidic KIA, H ⁺ and strongly basic AIA, Typ I, OH ⁻ 1:1,5
Kapazität (eq/l) capacity		min. 1.0
Trocknungsverlust (%) loss on drying	3h; 105°C	55.0 – 65.0
Korngröße (mesh) particle size		16 – 50 (0.3 – 1.2mm)
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.