
SERVA Tris-Tricin/SDS Proben Puffer (2x) SERVA Tris-Tricine/SDS sample buffer (2x)	Kat.Nr./Cat.No. :	42551
--	--------------------------	--------------

Parameter parameter	Spezifikation specification
Aussehen appearance	dunkelblau violette Lösung dark blue violet solution
pH	8.0 – 9.0
Zusammensetzung composition	
TRIS/ HCl pH 8.45	900mM
Glycerin/ glycerol (#23176)	24%
SDS (# 20763)	4%
SERVA Blau G/ Blue G (#35050)	0.015%
Phenolrot/ phenol red (#32095)	0.005%
Anwendung application	für SDS-PAGE for SDS-PAGE
Lagerung (°C) storage	+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.