

Collagen R Solution 0.4%
sterile

Cat.No.: 47256

Parameter	Method	Specification
Molecular weight		ca. 300 000
Source		rat tail
Concentration		4 mg/ ml in 0.1% acetic acid
Sterility		passes test
Gelling test		corresponds to reference
Storage (°C)		+2 to +8
References		Strom, S. C. & Michalopoulos, G. (1982) Methods Enzymol. 82 , 544 - 55 Miller, E. J. (1976) Mol. Cell. Biochem. 13 , 165 – 91

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.