

Albumin bovine cell culture grade, lyophil.	Cat.No. : 47330
---	------------------------

Parameter	Method	Specification
Appearance	Visual inspection	white to yellow with tan to green cast lyophilisate
Solubility	4 % in H ₂ O	clear to slightly hazy, faint yellow-green solution
Protein (%)	Nitrogen analysis	min. 96
Purity (%)	Agarose zone electrophoresis	min. 96
pH	10 % in H ₂ O	6.5 – 7.5
Endotoxin (EU/mg)	Limulus Amebocyte Lysate	max. 3
Moisture (%)	Loss on drying	max. 5.0
APC, average (cfu/g)	Aerobic plate count	max. 1000
Storage (°C)		-15 to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.