

Amphotericin B Lösung sterilfiltriert Amphotericin B solution sterile filtered	Kat.Nr./Cat.No. :	47982
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₄₇ H ₇₃ NO ₁₇
Molekulargewicht molecular weight		924.1
Aussehen appearance		klare gelbe Lösung clear yellow solution
Konzentration (mg/ml) concentration	in H ₂ O	0.25
pH		7.0 – 9.0
Osmolarität (mOsmol/l) osmolarity		max. 20
Sterilttest sterility		negativ negative
A 1cm/ 326 nm	1:20 verdünnt / diluted in H ₂ O	min. 0.7
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.