

<b>Tylosin-tartrat Lösung (100x)</b> sterilfiltriert	<b>Kat.Nr./Cat.No. :</b>	<b>47988</b>
Tylosine-tartrate solution (100x) sterile filtered		

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		C <sub>45</sub> H <sub>77</sub> NO <sub>17</sub> ·C <sub>4</sub> H <sub>6</sub> O <sub>6</sub>
<b>Molekulargewicht</b> molecular weight		1054.2
<b>Aussehen</b> appearance		klare farblose Lösung clear colorless solution
<b>Konzentration (mg/ml)</b> concentration (mg/ml)	in PBS	5
<b>pH</b>		7.0 – 7.5
<b>Osmolarität (mOsmol/l)</b> osmolarity		max. 320
<b>Sterilttest</b> sterility		negativ negative
<b>A 1cm/ 292 nm</b>	1:100 verdünnt / diluted in H <sub>2</sub> O	> 1.0
<b>Lagerung (°C)</b> storage		-15 bis/ to -25

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> We do not guarantee that the product can be used for a special application.</p>
--