

Pepton aus Casein tryptisch-leichtlöslich-für Gewebekultur Peptone from casein, tryptic, readily soluble	Kat.Nr./Cat.No. :	48600
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		hellelfenbeinfarbenes Pulver light-ivory colored powder
Trocknungsverlust (%) loss on drying		max. 6.0
Asche (%) ash		max. 15.0
pH	2 % in H ₂ O	6.5 – 7.5
Gesamtstickstoff (%) total nitrogen		min. 10
Aminostickstoff (%) amino nitrogen		min. 3.9
Mikrobiologische Analyse microbiological analysis	Salmonellen salmonella	negativ negative
	stand. plate count (col/ g)	< 5 000
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.