

Pepton aus Fleisch, pancreatic Peptone from meat, pancreatic	Kat.Nr./Cat.No. :	48619
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		beiges Pulver beige powder
pH	5 % in H ₂ O	5.0 – 6.0
Trocknungsverlust (%) loss on drying		max. 6.0
Gesamtstickstoff (TN) (%) total nitrogen		15.0 – 16.0
Aminostickstoff (AN) (%) amino nitrogen		3.0 – 4.0
AN/TN x 100		18 – 26
Sulfatasche (%) residue on ignition		max. 7.0
Chlorid (%) chloride	als/ as NaCl	max. 3.0
Mikrobiologische Analyse microbiological analysis	total aerobic microbial count	< 10 000/ 1 g
	Staphylococcus aureus	absence/ 10 g
	Salmonella	absence/ 25 g
	Escherichia coli	absence/ 1 g
	Coliformes	< 10/ 1 g
	Yeasts and moulds	< 20/ 1 g
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.