

---

<b>Pepton aus Gelatine</b> , pankreatisch Peptone from gelatin, pancreatic	<b>Kat.Nr./Cat.No. :</b>	<b>48622</b>
---	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Aussehen</b> appearance		cremefarbenes Pulver cream-colored powder
<b>Löslichkeit</b> solubility	5% in H <sub>2</sub> O	klare gelbliche Lösung clear yellowish solution
<b>Wasser (%)</b> water	Trocknungsverlust loss on drying	max. 6.0
<b>Asche (%)</b> ash		max. 8.0
<b>pH</b>	5% in H <sub>2</sub> O	6.5 – 7.5
<b>Gesamtstickstoff (TN) (%)</b> total nitrogen		15.5 – 16.5
<b>Aminostickstoff (AN) (%)</b> amino nitrogen		2.4 – 3.4
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.