

Pepton aus Lactalbumin pancreatisch, leichtlöslich Peptone from lactalbumin pancreatic, readily soluble	Kat.Nr./Cat.No. :	48625
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		hell elfenbeinfarbenes Pulver light-ivory colored powder
Trocknungsverlust (%) loss on drying		max. 6.0
Asche (%) ash		max. 15.0
pH	6 % in H ₂ O	6.7 – 7.7
Gesamtstickstoff (TN) (%) total nitrogen		min. 10.0
Aminostickstoff (AN) (%) amino nitrogen		min. 4.8
Mikrobiologische Analyse microbiological analysis		
	Salmonellen salmonella	negativ negative
	coliforms	negativ negative
	stand. plate count	< 5000/ g
Lagerung (°C) storage		+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--