

Trypton aus Casein, pankreatisch Tryptone from casein, pancreatic	Kat.Nr./Cat.No. :	48647
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		cremefarbenes Pulver cream colored powder
pH	2 % in H ₂ O	6.5 – 7.5
Trocknungsverlust (%) loss on drying		max. 6.0
Gesamtstickstoff (TN) (%) total nitrogen		12.5 – 13.5
Aminostickstoff (AN) (%) amino nitrogen		3.5 – 4.5
AN/TN x 100		26 – 37
Sulfatasche (%) residue on ignition		max. 15.0
Chlorid (%) chloride	als/ as NaCl	max. 1.0
Mikrobiologische Analyse microbiological analysis	total aerobic microbial count	< 10 000/g
	Staphylococcus aureus	absence /10g
	Salmonella	absence /25g
	Escherichia coli	absence /10g
	Coliforms	< 10/g
	Yeasts and moulds	< 20/g
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.