
Pepton aus Lactalbumin, pankreatisch hydrolysiert Peptone from lactalbumin, pancreatic hydrolysate	Kat.Nr./Cat.No. :	48662
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		grauweißes bis hellgelbes Pulver off-white to light yellow powder
Trocknungsverlust (%) loss on drying		max. 5.0
Gesamtstickstoff (TN) (%) total nitrogen		min. 12.0
Sulfatasche (%) residue on ignition		max. 7.5
pH	1 % in H ₂ O (25°C)	5.9 – 6.9
Chlorid (%) chloride	als/ as NaCl	max. 5.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.