

Pepton IV, pflanzlich Peptone from vegetable IV	Kat.Nr./Cat.No. :	48669
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		hellgelbes bis gelbes grünstichiges Pulver light yellow to yellow slight greenish powder
Trocknungsverlust (%) loss on drying		max. 4.5
Gesamtstickstoff (TN) (%) total nitrogen		min. 11.5
Aminostickstoff (AN) (%) amino nitrogen		min. 3.8
Sulfatasche (%) residue on ignition		max. 11.0
pH	1 % in H ₂ O (25°C)	6.3 – 7.3
Chlorid (%) chloride	als/ as NaCl	max. 4.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.