
Anti-Collagenase (<i>Clostridium histolyticum</i>) polyclonal antibody from sheep, solution	Cat. No. :	58050
---	-------------------	--------------

Parameter	Method	Specification
Purity	SDS-PAGE	> 85%
IgG concentration	OD (280nm)	1.0 mg/ml

For customer information:

Origin	IgG from sheep IgG purified by protein G affinity chromatography
Antigen	Collagenase from <i>Clostridium histolyticum</i>
Buffer	100 µl PBS, 0.01 % sodium azide
Storage conditions	-15 to -25 °C

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.