

L-Selenomethionin L-Seleno-Methionine	Kat.Nr./Cat.No. :	77765
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₅ H ₁₁ NO ₂ Se
Molekulargewicht molecular weight		196.1
Aussehen appearance		weißes, kristallines Pulver white, crystalline powder
Gehalt (%) assay		min. 99.0
Identität identity	IR	entspricht corresponds
Selen (%) selenium	AAS	39.0 – 41.0
Schwermetalle (ppm) heavy metals		max. 20
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.