

Uranylacetat ·2H₂O reinst Uranyl acetate·2H ₂ O research grade	Kat.Nr./Cat.No. :	77870
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		(CH ₃ COO) ₂ UO ₂ ·2H ₂ O
Molekulargewicht molecular weight		424.2
Aussehen appearance		gelbes Pulver yellow powder
Gehalt (%) assay	Titration	min. 98.0
Unlösliches (%) insolubles		max. 0.01
Chlorid (%) chloride	Cl ⁻	max. 0.003
Sulfat (%) sulfate	SO ₄ ²⁻	max. 0.01
Schwermetalle (%) heavy metals	als/ as Pb	max. 0.002
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.