

Dextrans

Dextrans fractionated, research grade

Dextrans of narrow molecular weight distribution.

Loss on drying max. 7 %

Sulfated ash max. 0.3 %

Heavy metals (Pb) max. 10 ppm

Type	Mol.-Weight	Cat.No.	Quantity	Cat.No.	Quantity
Dextran FP40	35 000 – 45 000	18665.01	100 g	18665.02	500 g

Dextrans fractionated, technical grade

Type	Mol.-Weight	Cat.No.	Quantity	Cat.No.	Quantity
Dextran 4	3 000 – 4 500	18687.01	100 g	18687.02	500 g
Dextran 8	8 000 – 12 000	18689.01	100 g	18689.02	500 g
Dextran 100	100 000 – 200 000	18693.01	100 g	18693.02	500 g
Dextran 500	ca. 500 000	18696.01	100 g	18696.02	500 g

Modified Dextrans

Type	Mol.-Weight	Cat.No.	Quantity
DEAE-Dextran research grade	Dextran: 150 000	18686.01	25 g
		18686.02	100 g
Dextran sulfate 500 - Na-salt stabilized (0.5 – 2 % phosphate)	500 000	18706.01	25 g
		18706.02	50 g
		18706.03	100 g

Ver 07/12