

Biochemicals
Electrophoresis
Bioseparation
Life Sciences
Specials

Hoefer PR250 Mini Orbital Shaker

Adjustable speed, compact and quiet

New: Mini Orbital Shaker

The new Hoefer PR250 Mini Orbital Shaker with adjustable speed is a compact and quiet instrument.

The digital LED display and easy programming sequence makes setting speed and time simple; the timer can be set from 0 to 99 hours 59 minutes or set to run continuously. The PR250 Mini Orbital Shaker speed can be adjusted to from 50 to 250 rpm, and has an orbit of 20 mm. The unit has a small footprint and is cold room compatible for applications that require refrigeration.

- Variable speed from 50 to 250 rpm
- Cold room compatible
- Flat platform with skid resistant mat
- Orbit 20 mm
- Quiet sealed bearing rotation
- LED display
- 99 hours 59 minutes timer
- Audible end of run alert
- Capable of continuous run
- Small footprint
- Optional accessory platforms available



PR250 Technical Specifications

Platform	23 x 27 cm
Speed	50 to 250 rpm
Timer	0 to 99 h 59 min, continuous
Unit Dimensions (w x h x d)	28 x 11 x 27 cm
Weight	2.5 kg
Universal Voltage Input	100 - 240 VAC, 50/60 Hz
Safety Certifications	EN61010-1:2001, CE

Each model includes: Mini Orbital Shaker, flat platform with skid resistant mat, screw driver.

There are three accessory tables available separately which can be easily installed: a twelve place 100 ml beaker clamping table, a six place 250 ml beaker clamping table and a universal table with two adjustable bar clamps.

Ordering Information

Cat. No.	Product
PR250	Mini Orbital Shaker
PR250-F100	Platform for 12 x 100 - 150 ml flasks
PR250-F100	Platform for 6 x 250 - 300 ml flasks
PR250-UTBL	Universal Platform with adjustable bars
PR250-FTBL	Flat Platform with skid resistant mat