

Biochemicals
Electrophoresis
Bioseparation
Life Sciences
Specials

SERVA IPG *BlueStrips*


New strip lengths and pH gradient

Broadest range of IPG strips now available from SERVA!


With the latest addition of two new strip lengths (11 and 17 cm) and one new pH gradient (5 - 8) SERVA offers with its SERVA IPG *BlueStrips* the broadest range of IPG strips for the proteomics research community.

The new SERVA IPG *BlueStrips*:

- Strip length available: 7, 11, 17, 18 and 24 cm
- pH gradient available: 3 - 10, 3 - 10 NL, 3 - 6, 4 - 7, 5 - 8, 6 - 10
- For rehydration in convenient 2 ml pack sizes: SERVALYT™ Carrier Ampholytes 3 - 10, 3 - 10 IsoDalt, 3 - 6, 4 - 7, 5 - 8 and 6 - 9
- All gel strips are 3 mm wide
- Other pH gradients and gel strip lengths are available upon request



Overview: IPG strip lengths and pH gradients available from SERVA



The schematic drawing shows all strip lengths and pH gradients available for SERVA IPG *BlueStrips*, the new pH gradient and strip lengths are indicated in red. All combinations of strip length and pH gradient are available, allowing the researcher to choose among 30 different types of IPG strips. For rehydration convenient 2 ml pack sizes of SERVALYT™ carrier ampholytes are available to be added to the rehydration buffer for optimum results in your demanding 2D gel electrophoresis experiments.

Product	SERVA IPG <i>BlueStrips</i> (12 strips/pack)					
pH Range	3 - 10	3 - 10 NL	3 - 6	4 - 7	5 - 8	6 - 10
Length 7 cm	43001.01	43002.01	43005.01	43003.01	43006.01	43004.01
Length 11 cm	43031.01	43032.01	43035.01	43033.01	43036.01	43034.01
Length 17 cm	43041.01	43042.01	43045.01	43043.01	43046.01	43044.01
Length 18 cm	43011.01	43012.01	43015.01	43013.01	43016.01	43014.01
Length 24 cm	43021.01	43022.01	43025.01	43023.01	43026.01	43024.01
Product	Recommended for rehydration: SERVALYT™ (2 ml/pack)					
pH Range	3 - 10	3 - 10 IsoDalt	3 - 6	4 - 7	5 - 8	6 - 9

