

Product of the Month

SERVA
■ serving scientists ■

Get **20% discount*** now!

March

*For all SERVA customers within the EU. Outside Germany and Austria online shop orders only. Can not be combined with other discounts, promotions etc.

Tergitol™ 15-S-9 – Biodegradable Alternative for Triton X-100

Tergitol™ 15-S-9

Tergitol™ 15-S-9 is a biodegradable alternative for Triton X-100. It is a non-ionic detergent containing a secondary ethoxylated alcohol resulting in a "T"-shaped molecular structure (Fig.1). The detergent is compatible with anionic, cationic, and other non-ionic detergents.

Applications

- Isolation, purification and analysis of membrane components, e.g. lipid rafts
- Solubilization of hydrophobic proteins, e.g. integral membrane proteins like phosphatidylinositol kinase




Fig. 1: General structure and formula of Tergitol™ 15-S-9

Features

- Rapid dissolution
- Compatible with anionic, cationic and non-ionic detergents
- Readily biodegradable (OECD 301E)
- Low foam stabilities
- Low odor



Tergitol = Trademark of DOW Chemical Company

Ordering Information

	Cat. no.	Size	List price	Special price*
Tergitol™ 15-S-9	37242.01	100 ml	12.00 €	9.60 €
	37242.02	500 ml	28.00 €	22.40 €
	37242.03	2.5 L	67.00 €	53.60 €

*** The special prices are valid until March 31st, 2021.**

For further information, please contact:

Dr. Judith Koch
Tel.: +49 6221 13840-44
E-Mail: Tech.Service@serva.de

SERVA Electrophoresis GmbH
Carl-Benz-Str. 7
D-69115 Heidelberg • Germany
www.serva.de